

Product datasheet for **GA116939**

Human SEC14L3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEC14L3 CRISPRa kit - CRISPR gene activation of human SEC14 like lipid binding 3
Symbol:	SEC14L3
Components:	GA116939G1, SEC14L3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116939G2, SEC14L3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116939G3, SEC14L3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UDX4</u>
Synonyms:	TAP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)