

## Product datasheet for **GA116950**

### Human RGL4 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | RGL4 CRISPRa kit - CRISPR gene activation of human ral guanine nucleotide dissociation stimulator like 4  |
| Symbol:       | RGL4  |
| Components:   | GA116950G1, RGL4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA116950G2, RGL4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA116950G3, RGL4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q8IZI4</a>  |
| Synonyms:     | Rgr   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)   |

#### Product images:



[View online »](#)