


Product datasheet for **GA116953**

Human GPR110 (ADGRF1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ADGRF1 CRISPRa kit - CRISPR gene activation of human adhesion G protein-coupled receptor F1
Symbol:	ADGRF1
Components:	GA116953G1, GPR110 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116953G2, GPR110 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116953G3, GPR110 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T601</u>
Synonyms:	GPR110; hGPCR36; KPG_012; PGR19
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)