

Product datasheet for **GA116955**

Human PGBD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PGBD3 CRISPRa kit - CRISPR gene activation of human piggyBac transposable element derived 3
Symbol:	PGBD3
UniProt ID:	Q8N328
Synonyms:	FLJ90201
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)