

Product datasheet for **GA116974**

Human **STK32C** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	STK32C CRISPRa kit - CRISPR gene activation of human serine/threonine kinase 32C
Symbol:	STK32C
UniProt ID:	Q86UX6
Synonyms:	PKE; YANK3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)