

Product datasheet for **GA117020**

Human MORN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MORN3 CRISPRa kit - CRISPR gene activation of human MORN repeat containing 3
Symbol:	MORN3
Components:	GA117020G1, MORN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117020G2, MORN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117020G3, MORN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PF18</u>
Synonyms:	MGC71497
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)