

Product datasheet for **GA117024**

Human **GAS2L3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	GAS2L3 CRISPRa kit - CRISPR gene activation of human growth arrest specific 2 like 3
Symbol:	GAS2L3
Components:	GA117024G1, GAS2L3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117024G2, GAS2L3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117024G3, GAS2L3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q86XJ1
Synonyms:	G2L3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)