

Product datasheet for **GA117029**

Human C12orf40 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C12orf40 CRISPRa kit – CRISPR gene activation of human chromosome 12 open reading frame 40
Symbol:	C12orf40
Components:	GA117029G1, C12orf40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA117029G2, C12orf40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA117029G3, C12orf40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q86WS4</u>
Synonyms:	HEL-206; HEL-S-94
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

