

Product datasheet for **GA117034**

Human SIAH3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SIAH3 CRISPRa kit - CRISPR gene activation of human siah E3 ubiquitin protein ligase family member 3
Symbol:	SIAH3
UniProt ID:	Q8IW03
Synonyms:	FLJ39203
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)