

## Product datasheet for **GA117035**

### Human **KCNRG** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNRG CRISPRa kit - CRISPR gene activation of human potassium channel regulator
Symbol:	KCNRG
UniProt ID:	<a href="#">Q8N5I3</a>
Synonyms:	DLTET
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)