


Product datasheet for **GA117038**

Human GPR137C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GPR137C CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor 137C
Symbol:	GPR137C
Components:	GA117038G1, GPR137C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117038G2, GPR137C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117038G3, GPR137C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N3F9</u>
Synonyms:	TM7SF1L2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)