

Product datasheet for **GA117044**

Human TSSK4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSSK4 CRISPRa kit - CRISPR gene activation of human testis specific serine kinase 4
Symbol:	TSSK4
UniProt ID:	Q6SA08
Synonyms:	C14orf20; STK22E; TSK-4; TSK4; TSSK5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)