

Product datasheet for **GA117078**

Human TSEN54 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSEN54 CRISPRa kit - CRISPR gene activation of human tRNA splicing endonuclease subunit 54
Symbol:	TSEN54
UniProt ID:	Q7Z6J9
Synonyms:	PCH2A; PCH4; PCH5; sen54; SEN54L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)