

Product datasheet for **GA117107**

Human Laminin (LAMA1) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | LAMA1 CRISPRa kit - CRISPR gene activation of human laminin subunit alpha 1 |
| Symbol: | LAMA1 |
| UniProt ID: | P25391 |
| Synonyms: | LAMA; PTBHS; S-LAM-alpha |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)