


Product datasheet for **GA117109**

Human C18orf26 (DYNAP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DYNAP CRISPRa kit - CRISPR gene activation of human dynactin associated protein
Symbol:	DYNAP
Components:	GA117109G1, C18orf26 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117109G2, C18orf26 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117109G3, C18orf26 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N1N2</u>
Synonyms:	C18orf26
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)