

## Product datasheet for **GA117111**

### Human **SIGLEC15** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	SIGLEC15 CRISPRa kit - CRISPR gene activation of human sialic acid binding Ig like lectin 15
Symbol:	SIGLEC15
UniProt ID:	<u><a href="#">Q6ZMC9</a></u>
Synonyms:	CD33L3; HsT1361; SIGLEC-15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)