

## Product datasheet for **GA117118**

### Human ZNF776 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | ZNF776 CRISPRa kit - CRISPR gene activation of human zinc finger protein 776                                                                                                                                                                                                                                  |
| Symbol:       | ZNF776                                                                                                                                                                                                                                                                                                        |
| Components:   | GA117118G1, ZNF776 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117118G2, ZNF776 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117118G3, ZNF776 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q68DI1</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | DKFZp686K10134; DKFZp781G1213; FLJ30886; FLJ38288                                                                                                                                                                                                                                                             |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)