

## Product datasheet for **GA117128**

### Human **LRRC68 (PPP1R37)** activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PPP1R37 CRISPRa kit - CRISPR gene activation of human protein phosphatase 1 regulatory subunit 37
Symbol:	PPP1R37
UniProt ID:	<u><a href="#">O75864</a></u>
Synonyms:	LRRC68
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)