

Product datasheet for **GA117149**

Human **SLC25A41** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A41 CRISPRa kit - CRISPR gene activation of human solute carrier family 25 member 41
Symbol:	SLC25A41
Components:	GA117149G1, SLC25A41 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117149G2, SLC25A41 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117149G3, SLC25A41 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N5S1</u>
Synonyms:	APC4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

