

Product datasheet for **GA117156**

Human ZNF875 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF875 CRISPRa kit - CRISPR gene activation of human zinc finger protein 875
Symbol:	ZNF875
Components:	GA117156G1, ZNF875 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117156G2, ZNF875 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117156G3, ZNF875 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P10072</u>
Synonyms:	ZNF875
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)