

Product datasheet for **GA117157**

Human TAF A3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF A3 CRISPRa kit - CRISPR gene activation of human TAF A chemokine like family member 3
Symbol:	TAF A3
UniProt ID:	Q7Z5A8
Synonyms:	TAF A-3; TAF A3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)