


Product datasheet for **GA117170**

Human ANKRD34A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ANKRD34A CRISPRa kit - CRISPR gene activation of human ankyrin repeat domain 34A
Symbol:	ANKRD34A
Components:	GA117170G1, ANKRD34A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117170G2, ANKRD34A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117170G3, ANKRD34A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q69YU3
Synonyms:	ANKRD34
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)