

Product datasheet for **GA117174**

Human ZNF326 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF326 CRISPRa kit – CRISPR gene activation of human zinc finger protein 326
Symbol:	ZNF326
Components:	GA117174G1, ZNF326 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117174G2, ZNF326 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117174G3, ZNF326 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5BKZ1</u>
Synonyms:	dJ871E2.1; ZAN75; Zfp326; ZIRD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

