

Product datasheet for **GA117180**

Human C20orf203 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C20orf203 CRISPRa kit - CRISPR gene activation of human chromosome 20 open reading frame 203
Symbol:	C20orf203
Components:	GA117180G1, C20orf203 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117180G2, C20orf203 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117180G3, C20orf203 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NBC4</u>
Synonyms:	RP11-503P21.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)