

Product datasheet for **GA117199**

Human DUSP28 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DUSP28 CRISPRa kit - CRISPR gene activation of human dual specificity phosphatase 28
Symbol:	DUSP28
Components:	GA117199G1, DUSP28 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117199G2, DUSP28 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117199G3, DUSP28 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4G0W2</u>
Synonyms:	DUSP26; VHP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

