

Product datasheet for **GA117218**

Human **RPUSD3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	RPUSD3 CRISPRa kit - CRISPR gene activation of human RNA pseudouridine synthase D3
Symbol:	RPUSD3
Components:	GA117218G1, RPUSD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117218G2, RPUSD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117218G3, RPUSD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P087</u>
Synonyms:	FLJ34707; FLJ37268; MGC29784
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)