

Product datasheet for **GA117219**

Human PRRT3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRRT3 CRISPRa kit - CRISPR gene activation of human proline rich transmembrane protein 3
Symbol:	PRRT3
Components:	GA117219G1, PRRT3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117219G2, PRRT3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117219G3, PRRT3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5FWE3</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)