

## Product datasheet for **GA117226**

### Human CRIPAK activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | CRIPAK CRISPRa kit - CRISPR gene activation of human cysteine rich PAK1 inhibitor                                                                                                                                                                                                                             |
| Symbol:       | CRIPAK                                                                                                                                                                                                                                                                                                        |
| Components:   | GA117226G1, CRIPAK gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117226G2, CRIPAK gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117226G3, CRIPAK gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | FLJ34443; MGC150528                                                                                                                                                                                                                                                                                           |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)