

Product datasheet for **GA117233**

Human RING finger protein 175 (RNF175) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RNF175 CRISPRa kit - CRISPR gene activation of human ring finger protein 175
Symbol:	RNF175
Components:	GA117233G1, RING finger protein 175 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117233G2, RING finger protein 175 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117233G3, RING finger protein 175 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N4F7</u>
Synonyms:	FLJ34190
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)