

Product datasheet for **GA117238**

Human SH3PXD2B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3PXD2B CRISPRa kit - CRISPR gene activation of human SH3 and PX domains 2B
Symbol:	SH3PXD2B
Components:	GA117238G1, SH3PXD2B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117238G2, SH3PXD2B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117238G3, SH3PXD2B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A1X283</u>
Synonyms:	FAD49; FTHS; HOFI; KIAA1295; TKS4; TSK4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)