

Product datasheet for **GA117246**

Human **SLC36A3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC36A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 36 member 3
Symbol:	SLC36A3
Components:	GA117246G1, SLC36A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117246G2, SLC36A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117246G3, SLC36A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q495N2</u>
Synonyms:	PAT3; TRAMD2; tramdorin2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)