

Product datasheet for **GA117247**

Human KIF4B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIF4B CRISPRa kit – CRISPR gene activation of human kinesin family member 4B
Symbol:	KIF4B
Components:	GA117247G1, KIF4B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117247G2, KIF4B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117247G3, KIF4B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q2VIQ3</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

