

Product datasheet for **GA117269**

Human ZNF789 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF789 CRISPRa kit - CRISPR gene activation of human zinc finger protein 789
Symbol:	ZNF789
Components:	GA117269G1, ZNF789 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117269G2, ZNF789 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117269G3, ZNF789 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5FWF6</u>
Synonyms:	zinc finger protein 789
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)