

Product datasheet for **GA117278**

Human TRIQK activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRIQK CRISPRa kit - CRISPR gene activation of human triple QxxK/R motif containing
Symbol:	TRIQK
UniProt ID:	<u>Q629K1</u>
Synonyms:	C8orf83; PRO0845; UPF0599
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)