

Product datasheet for **GA117281**

Human NKAIN3 activation kit by CRISPRa

Product data:

Product Type: CRISPRa Kits

Description: NKAIN3 CRISPRa kit - CRISPR gene activation of human sodium/potassium transporting ATPase interacting 3


Symbol: NKAIN3

Components: GA117281G1, NKAIN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)
GA117281G2, NKAIN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)
GA117281G3, NKAIN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)
1 CRISPRa-Enhancer vector, SKU GE100056
1 CRISPRa scramble vector, SKU GE100077

Synonyms: FAM77D

Format: 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)