

Product datasheet for **GA117290**

Human C9orf75 (TPRN) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TPRN CRISPRa kit - CRISPR gene activation of human taperin
Symbol:	TPRN
Components:	GA117290G1, C9orf75 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117290G2, C9orf75 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117290G3, C9orf75 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4KMQ1</u>
Synonyms:	C9orf75; DFNB79
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)