


Product datasheet for **GA117293**

Human C9orf150 (LURAP1L) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LURAP1L CRISPRa kit - CRISPR gene activation of human leucine rich adaptor protein 1 like
Symbol:	LURAP1L
Components:	GA117293G1, C9orf150 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117293G2, C9orf150 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117293G3, C9orf150 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IV03</u>
Synonyms:	bA3L8.2; C9orf150; HYST0841; LRAP35b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)