

Product datasheet for **GA117305**

Human P2Y8 (P2RY8) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | P2RY8 CRISPRa kit - CRISPR gene activation of human P2Y receptor family member 8 |
| Symbol: | P2RY8 |
| UniProt ID: | Q86VZ1 |
| Synonyms: | P2Y8 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)