

Product datasheet for **GA117313**

Human EIF4E3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF4E3 CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 4E family member 3
Symbol:	EIF4E3
UniProt ID:	Q8N5X7
Synonyms:	eIF-4E3; eIF4E-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)