

Product datasheet for **GA117316**

Human **VNIR2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	VNIR2 CRISPRa kit – CRISPR gene activation of human vomeronasal 1 receptor 2
Symbol:	VNIR2
Components:	GA117316G1, VNIR2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117316G2, VNIR2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117316G3, VNIR2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NFZ6</u>
Synonyms:	VIRL2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

