

Product datasheet for **GA117317**

Human **VNIR3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	VNIR3 CRISPRa kit – CRISPR gene activation of human vomeronasal 1 receptor 3 (gene/pseudogene)
Symbol:	VNIR3
Components:	GA117317G1, VNIR3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117317G2, VNIR3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117317G3, VNIR3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	VIRL3, FKSG46
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

