

Product datasheet for **GA117318**

Human **VNIR4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	VNIR4 CRISPRa kit – CRISPR gene activation of human vomeronasal 1 receptor 4
Symbol:	VNIR4
Components:	GA117318G1, VNIR4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117318G2, VNIR4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117318G3, VNIR4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z5H5</u>
Synonyms:	VIRL4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

