

Product datasheet for **GA117319**

Human **VNIR5** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	VNIR5 CRISPRa kit – CRISPR gene activation of human vomeronasal 1 receptor 5 (gene/pseudogene)
Symbol:	VNIR5
Components:	GA117319G1, VNIR5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117319G2, VNIR5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117319G3, VNIR5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z5H4</u>
Synonyms:	VIRL5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

