

Product datasheet for **GA117335**

Human SNAI3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SNAI3 CRISPRa kit - CRISPR gene activation of human snail family transcriptional repressor 3
Symbol:	SNAI3
Components:	GA117335G1, SNAI3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117335G2, SNAI3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117335G3, SNAI3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3KNW1</u>
Synonyms:	SMUC; SNAIL3; Zfp293; ZNF293
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)