

Product datasheet for **GA117336**

Human H3C15 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST2H3A CRISPRa kit - CRISPR gene activation of human histone cluster 2 H3 family member a
Symbol:	H3C15
UniProt ID:	Q71DI3
Synonyms:	H3/n; H3/o
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)