

## Product datasheet for **GA117356**

### Human KRTAP20-2 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                                        |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                           |
| <b>Description:</b>  | KRTAP20-2 CRISPRa kit – CRISPR gene activation of human keratin associated protein 20-2                                                                                                                                                                                                                                |
| <b>Symbol:</b>       | KRTAP20-2                                                                                                                                                                                                                                                                                                              |
| <b>Components:</b>   | GA117356G1, KRTAP20-2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117356G2, KRTAP20-2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117356G3, KRTAP20-2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <a href="#">Q3LI61</a>                                                                                                                                                                                                                                                                                                 |
| <b>Synonyms:</b>     | KAP20.2                                                                                                                                                                                                                                                                                                                |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                                |

#### Product images:

