

## Product datasheet for **GA117374**

### Human **KCNK18** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNK18 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 18
Symbol:	KCNK18
UniProt ID:	<a href="#">Q7Z418</a>
Synonyms:	K2p18.1; MGR13; TRESK; TRESK-2; TRESK2; TRIK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)