

Product datasheet for **GA117378**

Human CENATAC activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CCDC84 CRISPRa kit - CRISPR gene activation of human coiled-coil domain containing 84
Symbol:	CENATAC
Components:	GA117378G1, CENATAC gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117378G2, CENATAC gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117378G3, CENATAC gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86UT8</u>
Synonyms:	DLNB14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

