


Product datasheet for **GA117379**

Human TMEM225 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM225 CRISPRa kit - CRISPR gene activation of human transmembrane protein 225
Symbol:	TMEM225
Components:	GA117379G1, TMEM225 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117379G2, TMEM225 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117379G3, TMEM225 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6GV28</u>
Synonyms:	PMP22CD; PPP1R154
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)