

Product datasheet for **GA117399**

Human Polyserase3 (PRSS53) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRSS53 CRISPRa kit - CRISPR gene activation of human serine protease 53
Symbol:	PRSS53
UniProt ID:	Q2L4Q9
Synonyms:	POL3S; UNQ308
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)